



HEMATOLOGY

CELL COUNTER ANALYSERS — EVERY TEST VOLUME, COVERED

Right-size your hematology bench — from entry-level 3-part counters to premium high-throughput 5-part analyzers — every tier backed by genuine reagents, on-site installation and an engineer who already knows the machine.

Sysmex

Mindray

Horiba

Medonic

Peerless

Immunoshop

Ozelle

Qualcyte

Ozorate

HEAD OFFICE

G-43, Sumel Business Park - 6,
Dudheshwar Road, Shahibaugh,
Ahmedabad - 380004.

Ritesh Patel
9879208668

Dipen Parekh
9879479162

Pankaj Patel
9979535554

Ankit Patel
7567246266

HEMATOLOGY — 3-PART DIFFERENTIAL ANALYSERS

ENTRY-LEVEL · EVERYDAY CBC TESTING

Sysmex XQ-320

3-Part Differential Hematology Analyzer



- 20 parameters incl. dedicated neutrophil count — no compromise on clinical detail
- Just 16 µL sample needed — less blood, faster turnaround
- Up to 70 samples/hour keeps your walk-in volume moving

Sysmex XP-300

3-Part Differential Hematology Analyzer



- Compact footprint — fits any bench, any budget
- Simple one-touch CBC workflow your staff learns in a day
- The trusted Sysmex name patients and doctors already recognise

Medonic M32

3 Part Cell Counter | Haematology Analyzer



- Rock-solid everyday CBC accuracy, shift after shift
- Compact & easy to maintain — low total cost of ownership
- Backed by trusted Boule Medical engineering

AutoLab HA-360

3 Part Diff. Hematology Analyzer



- Dependable CBC testing for everyday lab volume
- Compact, easy-to-maintain design — less downtime
- Consistent, repeatable results your reports can stand on

Mission HA-360

Next Generation Hematology Analyzer



- Sheath flow impedance technology for cleaner counts
- Integrated hemoglobin measurement — one pass, full picture
- PLCC & PLCR give you an edge competitors' basic counters miss

LifeDx H3Pro

3 Part Hematology Analyzer with Unmatched Features



- CBC + 3-Diff, 21 parameters + 3 histograms — big-machine data, small-machine price
- Up to 70 samples/hour throughput
- Stores 5,00,000 results — never lose a patient history again

Ozelle EHBT-50

Compact Hematology Cell Counter



- Touch-panel operation with ID card scanning — fast patient onboarding
- Fast, reliable CBC results at the point of care
- Small footprint, big convenience for space-tight clinics

HEMATOLOGY — 5-PART DIFFERENTIAL ANALYSERS

MID-TIER · ADVANCED DIFFERENTIAL

Mindray BC-5150

5-Part Differential Hematology Analyzer



- Flow Cytometry + tri-Angle Laser Scatter + Chemical Dye Method — lab-grade differential accuracy
- Dedicated optical Basophil channel most 3-part counters can't offer

Mindray BC-5130


5-Part Differential Hematology Analyzer



- Space-conscious benchtop design — 5-part power, compact footprint
- Reliable everyday CBC testing your team can trust
- Low maintenance keeps service calls (and costs) down

Mindray BC-5000


5-Part Differential Hematology Analyzer



- Built for growing labs stepping up from 3-part counters
- Simple operation, low maintenance — an easy sell to any lab owner
- Consistent, accurate CBC results, every single run

Mindray BC-5140

5-Part Differential Hematology Analyzer



- Flow Cytometry + tri-Angle Laser Scatter + Chemical Dye Method
- Dedicated optical Basophil counting channel for confident abnormal flagging

Medonic M51

CBC + 5 Part Differential Cell Counter



- Engineered for uptime — the analyzer that doesn't let your lab down
- Low maintenance needs mean lower cost-per-test, higher margin for you

AutoLab HA-560

5 Part Diff. Hematology Analyzer



- 3D scattergram visualisation doctors trust at a glance
- 29 parameters including ALY#, ALY%, IG#, IG% — genuinely advanced diagnostics
- Intelligent flagging catches abnormal samples before they slip through

Horiba Microsemi CRP LC-667G

Hematology Analyzer with Integrated CRP




- CBC + CRP on one platform — one machine, one AMC, twice the revenue stream
- Ideal for point-of-care & small labs wanting to add inflammation markers
- Fast, reliable combined results in a single visit

HEMATOLOGY — PREMIUM & ESR

HIGH-THROUGHPUT · PREMIUM 5-PART SERIES

Mindray BC-700 Series

*High-Throughput 5-Part
Differential Analyzer*



- All-in-one CBC + ESR — one analyzer, two revenue lines
- Built for busy reference labs

Mindray BC-6800Plus

*Premium 5-Part Differential
Analyzer*



- Body fluid mode expands your test menu instantly
- Advanced flagging cuts manual smear reviews

Mindray BC-6200

*Premium 5-Part Differential
Analyzer*



- Built for the medium-to-high volume segment
- Flagging accuracy your pathologists will trust

Mindray BC-6000

*Premium 5-Part Differential
Analyzer*



- Compact footprint, premium performance
- Low maintenance — premium tech, no premium hassle

Horiba Yumizen H500 & H550

*Premium 5-Part Differential
Analyzer*



- Open Tube, Closed Tube & Auto Loader — scale as you grow
- Touchscreen interface your staff adopt fast

ESR ANALYSERS · COMPLEMENTARY ADD-ON

Qualcyte 10

Automated ESR Analyser



- 10-sample batch processing keeps mornings moving
- An easy add-on sale alongside any analyzer

Ozorate 40

Automated ESR Analyzer



- 40-sample capacity for higher-volume labs
- Reliable ESR results clients ask for by name