



Peerless




Accre 8

An Automatic
Chemiluminescence
Immunoassay
(CLIA) Platform



- Compact & easy to operate
- Ready to use single dose cartridge
- No extra consumables required

Genuine quality speaks for itself






EXCELLENT AS CLIA

- ALP-APS-5 Chemiluminescence
- Excellent reagent performance
- Nano magnetic particles for separation
- Broad menu with featured items
- Reliable calibration and quality control system

EASY AS POCT

- Single-dose cartridge, ready-to-use
- Avoid cross-contamination
- Short turn-around-time (TAT)
- All-in-one design, no consumables and extra cost
- Built-in scanner reads reagent information automatically

ASSAY PROCEDURE



ACCURATE
CONVENIENT
CONSISTENT
RAPID
ECONOMICAL

Accre 8

Automatic Chemiluminescence
Immunoassay (CLIA) System



Space-saving

- Mini size bench-top analyzer
- Only 42 kg of weight



Ease of use

- 15.6 inch touch screen
- Delicate ergonomic design



Compact machine design

- No liquid path inside
- Innovative integrated design of magnetic separation and pipetting system
- Simplified mechanical movement design

System	Bench-top CLIA analyzer
Methodology	Chemiluminescence Immunoassay (CLIA)
Dimension	500 mm x 370 mm x 600 mm
Weight	42 kg
Throughput	Upto 36 tests / hour
Result	10 mins for the first result
Mode	Batch mode
Operation	15.6 inch colorful LCD touch screen
Channel	6 channels
Sample Type	Whole blood / Capillary blood / Serum / Plasma
Maintenance	24 hours standby, 600 days MTBF [®] , no consumable needed for maintenance
Accessories	Built-in thermal printer and barcode scanner
Data Connection	LIS / HIS connection RS232, USB, Internet connection

®MTBF : Mean Time Between Failure

TEST PANEL



Thyroid

- TSH
- T3
- T4
- Ft3
- Ft4
- Anti-Tg
- Anti-TPO



Infectious Disease

- HBsAg
- Angi-HBs
- HBeAg
- Anti-HBe
- Anti-HBc
- Anti-HCV
- Anti-TP
- HIV Combo
- Influenza A & B
- Chagas
- H Pylori (HP)
- Dengue NS1 Ag
- Dengue IgG / IgM Ab



Anaemia

- VB12
- Folate
- Ferritin



Bone Metabolism

- 25-OH Vitamin D Total
- PTH
- β -CrossLaps
- Total P1NP
- N-MID Osteocalcin



Fertility

- AMH
- LH
- Prolactin
- FSH
- Progesterone
- Testosterone
- Estradiol (E2)
- Total β -HCG
- Estriol
- Inhibin B
- DHEA-S
- ACTH
- Cortisol
- Aldosterone
- Direct Renin
- hGH



Cardiac

- hs cTnI
- hs CTnT
- cTnI
- BNP
- NT-proBNP
- CK-MB
- MYO
- H-FABP
- D-Dimer
- Lp-PLA2
- ST2
- MPO
- Homocysteine
- Copeptin
- GDF-15
- Galectin-3



Tumor

- AFP
- CEA
- TPSA
- CA 125
- PG I
- PG II
- Gastrin-17 (G-17)
- FPSA
- CA 15.3
- CA 19.9
- S100
- ProGRP



TORCH

- Toxo IgG
- Toxo IgM
- Rubella IgG
- Rubella IgM
- CMV IgG
- CMV IgM



Inflammation

- PCT
- IL-6
- FR-CRP
- ASO
- SAA



Glyco Metabolism

- C-peptide
- Insulin



Allergy

- IgE

▪ Available

▪ In Development

▪ Coming Soon

Imported & Marketed by:



Peerless Biotech Pvt. Ltd.

36/37, 'Millennium House'
M. K. Srinivasan Nagar Main Road
Rajiv Gandhi Salai, Perungudi
Chennai - 600 096, Tamil Nadu, India

☎ +91 44 2496 0541
✉ info@peerlessbio.com
🌐 www.peerlessbio.com

Customer care ☎ +91 93812 10866

Genuine quality speaks for itself